

25
CELEBRATING 25 YEARS

1993-2018

AQUA-NOVA

Probably the best Life-Cycle Cost in the world, since 1993

Laboratory Water For Injection & Pure Steam



ADVANTAGES


- One control cabinet
- One PoU of media
- One kit of spare parts
- One installation

SAVE ON


- Investment
- Space
- Validation
- Heat emission



WFI



PS



Tailor Built

Ready for production at any capacity

Acceptance FAT testing

All units are thoroughly tested for production, consumption and function during the factory pre-delivery testing.

The distillate quality is tested by an independent laboratory. Typical quality values are:

Conductivity, <0.5 µS/cm
 Cfu/ml, <0.1
 EU/ml, <0.125

Typical Range

LABORATORY STILLS, electrically heated

E-4-C	3 kW	4 L/h
E-8-C	6 kW	8 L/h
E-25-C	18 kW	25 L/h
E-50-C	35 kW	50 L/h
E-100-3	6 kW	100 L/h